

4 Wheel Rollator | Technical Spec Sheet

Product Overview

- Description:** This product is a premium quality rollator designed to assist patients with improved mobility.
- Colours:** Available in blue colour only
- Handles:** Both handles are height adjustable
- Seat:** Fitted with a wide, comfortable seat and back rest
- Basket:** A secure detachable basket for storing belongings



Standard Size

Description

Structure	Aluminium
Handle	Aluminium
	Moulded PVC grip
Back Rest	Aluminium with foam cover
Seat	Comfortable, durable foam with zip cover

Measurements

Total Length	54cm	Seat Width	34cm
Total Width	60cm	Seat Depth	31cm
Total Height	82 - 97.5cm	Back Height	20cm
Front Wheels	19cm	Weight Limit	100kg
Back Wheels	19cm		

Wide Size

Description

Structure	Steel	*
Handle	Steel	*
	Moulded PVC grip	
Back Rest	Aluminium with foam cover	
Seat	Comfortable, durable foam with zip cover	

Measurements

Total Length	54cm	Seat Width	46cm	*
Total Width	70cm	Seat Depth	33cm	*
Total Height	75 - 91cm	Back Height	20cm	*
Front Wheels	20cm	Weight Limit	136kg	*
Back Wheels	20cm			*

* Different to the Standard Size