

Adjustable Reciprocal Folding Two Step Walker | Technical Spec Sheet

Product Overview

Description: An easy-to-use reciprocal folding walker with stand assist
Colours: Available in silver colour only
Handles: Soft grip handles



Specifications

Description

An easy-to-use reciprocal folding walker with stand assist. The reciprocal walking frame walks with the user. There is no need to lift the walker forward with each step, rather this walker, when used correctly, steps with the user. This walker may also be used as a rigid frame.

This walker can be folded with the push of a quick single release button allowing for efficient storage and travel. The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Dimensions

Length	53cm
Width	47cm
Height	73 - 85cm (Adjustable)
Foot	Fitted with 27mm foot plugs at lower end

Weight Limit

Weight Limit	110kg
--------------	-------

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, such as using appropriate footwear and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.