

Freely

Mobility Aids
Product Catalogue





Freely

Overview

Welcome to Freely and our range of mobility aids designed to enhance user independence and improve quality of life. The Freely range includes crutches, walking sticks, walkers and rollators to meet specific user needs and provide the utmost comfort and support.

At Freely we understand the significance of mobility assistance in maintaining an active and fulfilling lifestyle. Whether someone is recovering from an injury, living with a long-term condition, or simply seeking additional stability and confidence while on the move, our range of mobility aids has been carefully chosen to cater to their unique requirements.

From crutches that offer exceptional stability and weight distribution to walkers and rollators that combine durability with effortless manoeuvrability, each product in the range has been tried and tested to meet stringent safety standards and has been designed for the most important user requirements - ergonomics, lightweight materials, stability, anti-slip features, shock absorption, adjustability, portability and visual appeal.

A Trusted Mobility Partner



Trusted Across SADC

At Freely Mobility, we don't just deliver products, we support progress, independence, and long-term customer trust. Whether you're a healthcare provider, distributor, or stockist, partnering with Freely means aligning with a brand focused on consistent quality, customer comfort, and business reliability.



Crutches

Elbow Crutch

An ergonomically designed elbow crutch to assist users with recovery, improve mobility and provide support for lower extremity injuries. The Elbow Crutch is sold in pairs and extensions are available (55cm in length). A double adjustable option is also available, allowing both the shaft and the cuff height to be adjusted. The double adjustable crutch is subject to MOQ's.



**Adult
Standard Grip**

**Adult
Soft Grip**

Colour Silver (Natural) Anodised

Weight Limit

Adult 120kg
Youth 80kg
Paediatric 60kg

Foot size

Adult 21mm
Youth & Paediatric 19mm

Adjustable Height

Adult 70,3 - 92,7cm
Youth 51 - 73cm
Paediatric 41 - 63cm

**From handle to ground*

Permanent Crutch

The Permanent Crutch is a dependable solution for long-term walking support. The Permanent Crutch can be cut to size and is sold in pairs.



Colour Silver (Natural) Anodised

Weight Limit

Standard 120kg
Heavy Duty 150kg

Foot size 24mm

Height 120cm

**From cuff to ground*

Gutter Crutch

The Gutter Crutch is designed for individuals with arthritis, impaired grip, or partial weight-bearing conditions. This crutch features a padded forearm trough for even weight distribution providing stability for users. The adjustable hand grip ensures a comfortable and secure hold. This product is sold in single units.



Colour Silver (Natural) Anodised

Foot size 22mm

Adjustable Height

110-120cm

**From handle to ground*

Details

- Forearm padded with sponge and Velcro strap
- Adjustable handle from 10-20cm
- Height 100-120,4cm

**From handle to ground*

Crutches

Aluminium Axilla Crutches

The Aluminium Axilla Crutch is an underarm crutch designed for short term use as a result of a lower-limb injury and is available in five sizes. The Aluminium Axilla Crutch has soft rubber underarm padding and hand grips. Each size option is height adjustable and the handle position can also be adjusted for a more precise user fit. This product is sold in pairs.



Colour Silver (Natural) Anodised

Weight Limit

XS-S 80kg
M-L 100kg
XL 114kg

Adjustable Height

XSmall 84-105cm
Small 94-115cm
Medium 115-135cm
Large 135-155cm
XLarge 145-165cm

**From underarm to ground*

Wooden Axilla Crutches

The Wooden Axilla Crutch is an underarm crutch designed for short term use as a result of a lower-limb injury and is available in three sizes. The Wooden Axilla Crutch has soft rubber underarm padding and hand grips. Each size option is height adjustable and the handle position can also be adjusted for a more precise user fit. This product is sold in pairs.



Colour Light Wood

Weight Limit

S 80kg
M 100kg
L 100kg

Adjustable Height

Small 96-117cm
Medium 110-130cm
Large 126-146cm

**From underarm to ground*

Walking Sticks

Classic Walking Stick

The perfect solution for the elderly, the physically challenged, or people who have trouble walking or standing. The Classic Walking Stick design is suitable for left and right handed users, and is sold in single units.



Colour

Bronze Anodised
Silver (Natural) Anodised
White Powder Coated
Black Powder Coated

Weight Limit 110kg

Foot size 19mm

Adjustable Height 70 - 90cm

**From handle to ground*



Quadruped Walking Stick

The perfect solution for the elderly, the physically challenged, or people who have trouble walking or standing. The design is suitable for left and right handed users, and is sold in single units.



Colour Silver (Natural) Anodised

Weight Limit 110kg

Foot size 16mm

Adjustable Height 70 - 90cm

**From handle to ground*

Difference between the Classic and Quadruped Walking Stick

While the Classic Walking Stick relies on a single point of contact, the Quadruped Walking Stick offers four points of contact, providing increased stability and balance. The wider base of the Quadruped Walking Stick offers more support and helps distribute weight evenly, reducing the risk of falls or loss of balance..

Walking Sticks

Folding Walking Stick

The perfect solution for the elderly, those with mobility challenges, or difficulty walking. Suitable for both left and right-handed users, its foldable design adds convenience, making it easy to take anywhere, and is sold in single units.



Colour Black Anodised

Weight Limit 110kg

Folded Length 30cm

Foot size 19mm

Adjustable Height 84-94cm

**From handle to ground*



S-Bend Walking Stick

The S-bend Walking Stick is a versatile aid designed to support walking and assist in maintaining good posture. This product is sold in single units.



Colour Silver (Natural) Anodised

Weight Limit 100kg

Foot size 22mm

Adjustable Height 72-94cm

Walking Pole

Designed for individuals transitioning from a walker, cane or crutch to a more independent gait, such as patients rehabilitating from injury or conditions like stroke. This product is sold in single units. Additional feet can be bought separately.



Colour Silver (Natural) Anodised

Weight Limit 120kg

Foot size 24mm

Standard Height 152cm

**Length can be cut to size*

Folding Walkers & Rollators

Folding Walker/Rollator Two Push Button

An easy-to-use folding walker designed as an ideal and dependable solution for those who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick double release button allowing for efficient storage and travel. The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.



Colour Silver (Natural) Anodised & Red Powder Coated

Weight Limit 110kg

Foot size 27mm

Adjustable Height

Adult 80 - 87cm

Youth 66 - 75cm

Paediatric 56 - 63cm

**From handle to ground*



Wheel Diameter 12,5cm



Folding Rollator

Wheels can be added to convert the walker into a rollator designed to be used primarily by people who are recuperating from leg or back injuries as well as hip and knee replacements. The two wheels at the front make for easier and more energy efficient propulsion whilst still providing a stable base of support.

Folding Walkers & Rollators

Reciprocal Folding Two Step Walker

An easy-to-use, adjustable Reciprocal Folding Walker with stand assist. The reciprocal walking frame walks with the user. There is no need to lift the walker forward with each step, rather this walker, when used correctly, steps with the user. This walker may also be used as a rigid frame.

This walker can be folded with the push of a quick single release button allowing for efficient storage and travel. The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Dimensions

Length 53cm, Width 47cm, Height 73 - 85cm



Colour Silver (Natural) Anodised & Black/Red Powder Coated

Weight Limit 110kg

Foot size 27mm

Folding Walker One Push Button

An easy-to-use Folding Walker designed as an ideal and dependable solution for those who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick single release button allowing for efficient storage and travel.

The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Dimensions

Length 58cm, Width 49cm, Height 79 - 97cm



Colour Silver (Natural) & Grey Anodised

Weight Limit 110kg

Foot size 27mm

Four Wheel Folding Rollator

A quality rollator designed to assist users with improved mobility. With an ergonomic design the lightweight frame is fitted with four large, sturdy castor wheels for efficient movement in every direction. Secure light touch hand brakes on both handles make stopping effortless, and both can be locked in place to ensure safety in a seated or standing position. The rollator folds flat to a compact size for easy storage or transport.

Details

- Both handles are height adjustable
- Wide, comfortable seat and back rest (height 56cm)
- Standard seat (31.5 x 34.5cm), Wide seat (33.5 x 46.5cm)
- Detachable basket for storing belongings
- 18cm diameter wheels for both standard and wide

Colour Silver (Natural) & Burgundy Anodised

Weight Limit

Standard 100kg

Wide 136kg



Folding Walkers

Folding Walker (Youth)

Easy-to-use and designed specifically for children who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick single release button allowing for efficient storage and travel.

The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Colour Yellow & Black Powder Coated

Dimensions

Length 47cm, Width 58cm, Height 58 - 68cm

Foot size 27mm



Folding Walker (Kiddies)

Easy-to-use and designed specifically for children who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick single release button allowing for efficient storage and travel.

The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Colour Silver (Natural) Anodised & Red Powder Coated

Dimensions

Length 35cm, Width 39cm, Height 49 - 54cm

Foot size 27mm



Folding Knee Scooter with Rear Suspension

Provides stable, pain-free mobility for recovery from foot, ankle, or lower-leg surgery, or below-knee amputation. Its four wheels and rear shock absorber ensure a smooth ride indoors and outdoors.

A rear spring suspension adds comfort, while dual handbrakes offer safe control and parking. Quick-release levers adjust the handle height and fold the scooter flat for easy transport, with a detachable bag to keep essentials close.

Colour Blue, White & Black Powder Coated

Dimensions

- Length 85cm, Width 43cm, Height 84 - 99cm
- Knee Pad Length 32cm, Width 18cm, Height 49cm - 56cm

Unit weight 10kg

Wheel Diameter 20cm



Bathroom Mobility Aids

Shower Chair

Designed for patients that require extra stability and support for safer bathing. It features a lightweight aluminium frame, height-adjustable legs, a removable backrest, and padded armrests for easier side transfers. Non-slip feet keep it secure, and the drain-through seat reduces slipping.

Colour Silver (Natural) Anodised

Weight Limit 100kg

Dimensions

- Width 52cm, Depth 36cm, Height 77cm - 87cm
- Seat Width 40cm, Depth 33cm

Foot size 27mm



3 in 1 Basic Commode

Designed to work as a bedside commode, raised toilet seat, or sturdy safety frame. Its durable powder-coated steel frame offers reliable support, while moulded armrests and a fixed backrest add comfort and stability. Tool-free height adjustment and a 5L slide-out bucket with lid make cleaning and daily use simple.

Colour White Powder Coated

Weight Limit 100kg

Dimensions

- Width 56cm, Height 64cm - 80cm
- Width between arms 49cm
- Seat Width 36cm, Depth 44cm
- Armrest height 27cm

Foot size 27mm



Brace & Support Products

Pneumatic Walker Boot

Designed to provide superior stabilisation, comfort, and protection for individuals recovering from left or right foot or ankle injuries. Featuring an advanced pneumatic (air) system, this boot offers customisable compression to aid healing, reduce swelling, and enhance patient comfort.

Colour Black with Grey sole

Sizes

Small
Medium
Large
XLarge



Ankle Brace

Designed to deliver support, stability, and comfort for active individuals and those recovering from injury. This versatile brace combines durability with lightweight performance, ensuring users stay protected without compromising mobility.

Colour Black

Sizes

XSmall
Small
Medium
Large
XLarge
2XLarge



Hinged Knee Wrap

Designed to provide targeted support, stability, and comfort for individuals recovering from left or right knee injuries or managing ongoing knee conditions. Ideal for active individuals, athletes, and everyday users, this brace combines durable materials with an ergonomic design to help protect and support the knee without restricting natural movement.

Colour Black

Sizes

XSmall
Small
Medium
Large
XLarge
2XLarge
3XLarge



Brace & Support Products

Arm Sling

Provides secure, comfortable, and adjustable support for the left or right arm, elbow, and wrist during recovery from fractures, sprains, or post-surgical procedures. Lightweight yet durable, it offers excellent immobilisation while ensuring all-day comfort for the patient.

Colour Black

Sizes

Small
Medium
Large
XLarge
2XLarge



Thumb Spica

Provides targeted support and immobilisation for the wrist and thumb. Crafted from high-quality breathable materials, with integrated aluminium stays for stabilisation, it ensures comfort while effectively aiding recovery from injuries and chronic conditions. Suitable for both left and right hands, it combines versatility, comfort, and durability in one universal solution.

Colour Black

Sizes

XSmall - Small
Medium - Large
XLarge - 2XLarge



Wrist Brace

Designed to provide effective support and relief for sprained, injured, or painful wrists, whether during daily activities, rest, or sleep. Manufactured to the highest quality standards, this brace combines comfort, durability, and adjustability, making it suitable for a wide range of users.

Colour Black

Sizes

XSmall - Small
Medium - Large
XLarge - 2XLarge



Paediatric Arm Sling

Specially designed for kids, this sling provides secure, comfortable, and adjustable support for the left or right arm, elbow, and wrist during recovery from fractures, sprains, or after surgery. The lightweight yet durable design ensures proper immobilisation, while the playful animal print adds a touch of fun and helps little patients feel more at ease during healing.

Mobility Accessories

Toilet Seat Raiser

Used by patients especially after hip replacement surgery or other situations requiring an elevated seated position.

Colour White

Weight Limit 140kg

Details

- Fits all standard size toilets
- Made from white plastic (LLDPE)
- Features stainless steel feet with rubber tips
- Non-adjustable



Standard Grey Cuffs

Designed for the standard adult elbow crutch these replacement cuffs are sold in pairs.

Colour Grey



Standard Feet

Standard feet made from black and grey PVC, a thick base and galvanised washer inside. This product is sold in pairs.

Details

- Should be replaced regularly
- Easy to attach
- Available in two colour options (specific sizes per colour)



Colour Grey

Foot sizes 19, 21, 27mm



Colour Black

Foot sizes 16, 22, 24mm

A Partner You Can Depend On

Why Distributors, Pharmacies, Hospitals and Clinics Choose Freely



Consistent Stock Availability

Locally made and distributed with reliable turnaround times.



Competitive Pricing

We help you stay profitable without compromising quality.



Bulk Ordering & Account Support

Work directly with a supplier who understands your needs.



Dependable Supply Chain

Streamline logistics and support to minimise friction.

Partner with Us

At Freely Mobility, we take pride in being a trusted supplier of mobility aids. We understand the importance of providing you with quality products at competitive prices that will positively impact your business. Our goal is to help you build a reputation for delivering reliable mobility aid products to your own customers, enhancing customer satisfaction, and overall competitiveness in the market.

Here are some key advantages of partnering with us:

Strong Partnerships

We prioritise building long-term partnerships fostering open communication, trust, and mutual growth.

Expertise and Experience

We have years of experience sourcing, manufacturing and supplying quality products throughout Southern Africa.

Wide Product Range

Our products cater to various end user needs, providing you with flexibility in your offering and we're constantly developing or sourcing new, innovative products.

High-Quality Products

We provide reliable and durable products ensuring consistent product quality for your customers.

Freely

 www.freelymobility.co.za

 info@freelymobility.co.za