

Product Overview

Description:	Elbow crutch comfortable for medium to long term use
Type:	Underarm, adjustable
Tube finish:	Anodised matt 15 microns (μ)
Colours:	Available in natural colour only
Plastics:	Grey
Structure:	Adjustable aluminium tube
Leg adjustment:	Adjustable to 9 different lengths (25mm apart)
Punch options:	Single; or double by arrangement

Classification: ISO9999

Code: 120306



General Information

Support Angle	90°
Grip Angle	90°

Fittings

Under Arm	Made of polypropylene
	Fitted with moulded rubber padding
Handle	Made of polypropylene
	Fitted with moulded rubber / PU foam
Foot	Design for maximum adherence
	Replaceable

Upper Support (Aluminium 6063-T6)

Outside Diameter	19.43mm
Height Adjustment	X 9 holes
Plastics	Fitted with plastic slide at lower end for noise and vibration reduction

Lower Tube (Aluminium 6063-T6)

Outside Diameter	21.6mm
Height Adjustment	Fitted with plated spring clip
Plastics	Fitted with foot at lower end
Foot	Fitted with 22mm foot plug at lower end

Key Differences

	X-Small	Small	Medium	Large	X-Large
Mass	649g	725g	793g	877g	904g
Max Weight of User	80kg	80kg	100kg	100kg	114kg
Leg Section Length	84 - 105cm	94 - 115cm	115 - 135cm	135 - 155cm	145 - 165cm
Forearm Support Length	750mm	830mm	90°	90°	90°
Lower Tube Foot Length	345mm	345mm	345mm	345mm	440mm
Lower Tube Mass	145g	145g	145g	145g	173g

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, such as using appropriate footwear and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.

Please get in touch for more information

www.freelymobility.co.za @ info@freelymobility.co.za