

Ankle Brace | Technical Spec Sheet

Product Overview

Description:	Adult Ankle Brace
Type:	Adjustable
Colours:	Black
Materials:	Nylon-polyester mesh body with wide elastic nylon straps and hook-and-loop (Velcro) closures, breathable foam lining.



Description

Designed to deliver support, stability, and comfort for active individuals and those recovering from injury. This versatile brace combines durability with lightweight performance, ensuring users stay protected without compromising mobility.

Sizes

Small	30 - 33cm
Medium	33 - 36cm
Large	36 - 39cm
XLarge	39 - 41cm
2XLarge	41 - 45cm

*Size measured by the circumference around the ankle, taken from the base of the heel over the malleoli (ankle bone).

Indications

- Post-injury rehabilitation and immobilisation
- Support for sore or weak ankles
- Prevention or treatment of acute ankle sprains
- Ligament and tendon injuries

Features

- Quick-lace system for easy application and removal – more convenient than traditional lace-up braces.
- High-quality, durable materials resistant to pilling and wear.
- Secure fit for reliable ankle stabilisation.
- Ambidextrous design—fits left or right ankle.

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.