

Paediatric Arm Sling | Technical Spec Sheet

Product Overview

Description:	Paediatric arm sling
Type:	Adjustable
Colours:	Blue with pattern design
Materials:	Breathable, soft material construction with hook-and-loop (Velcro) closures for adjustable support and comfort.



Description

Specially designed for kids, this sling provides secure, comfortable, and adjustable support for the left or right arm, elbow, and wrist during recovery from fractures, sprains, or after surgery. The lightweight yet durable design ensures proper immobilisation, while the playful animal print adds a touch of fun and helps little patients feel more at ease during healing.

Indications

- Stabilises the shoulder, arm and wrist after injury or surgery.
- For fractures, rotator cuff injuries, sprains and dislocations.

Features

- Breathable design for airflow and reduced heat build-up.
- Adjustable shoulder strap for a quick, personalised fit.
- Removable swathe strap for added immobilisation.
- Soft foam neck/shoulder strap for enhanced comfort.
- Ambidextrous design—fits left or right arm.
- Hook-and-loop closures for easy on/off.
- Lightweight design for day-long wear.
- Machine washable to keep clean and fresh.

Sizes

	Sling Length	Max Sling Depth	Neck Strap Length	Waist Strap Length
One Size	23.5cm	13.5cm	40cm	40cm

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.

Please get in touch for more information

www.freelymobility.co.za @ info@freelymobility.co.za

Ver 1 (5 Dec 2025)