

Thumb Spica | Technical Spec Sheet

Product Overview

Description:	Adult Thumb Spica
Type:	Adjustable
Colours:	Black
Materials:	Breathable neoprene, nylon and polyester padded construction with laced as well as hook-and-loop (Velcro) closures.



Description

Provides targeted support and immobilisation for the wrist and thumb. Crafted from high-quality breathable materials, with integrated aluminium stays for stabilisation, it ensures comfort while effectively aiding recovery from injuries and chronic conditions. Suitable for both left and right hands, it combines versatility, comfort, and durability in one universal solution.

Sizes

	Palm	Forearm
XSmall - Small	17 - 21cm	20 - 24cm
Medium - Large	21 - 25cm	24 - 28cm
XLarge - 2XLarge	25 - 29cm	29 - 34cm

Indications

- Fractures, sprains and strains
- Thumb sprains and strains – especially of the ulnar collateral ligament (e.g. “Skier’s thumb” or “Gamekeeper’s thumb”).
- Arthritis
- Post-surgical or post-cast support

Features

- Breathable design for airflow and reduced heat build-up.
- Adjustable (Velcro) straps for a quick, personalised fit.
- Soft foam wrist and thumb padding for enhanced comfort.
- Lightweight design for day-long wear.
- Suitable for day- and night-time wear.
- Thumb loop for easy application and removal; hook-and-loop strap for quick one-handed fastening.
- Ambidextrous design—fits left or right hand.
- Durable build for long-term use.

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.