

Bariatric Folding Walker | Technical Spec Sheet

Product Overview

Description: An easy-to-use folding walker
Colours: Available in silver colour only
Handles: Soft grip handles



Specifications

Description

Provides a practical, dependable solution for those requiring stable support during daily activities. To ensure an ergonomic fit, the walker features 8 adjustable height settings via a simple one-push-button mechanism on each leg. Soft-grip handles are included to enhance comfort and reduce hand pressure during extended use.

For added convenience, a two-push-button release system allows the unit to fold quickly for storage or travel. While designed as a standard folding walker, it is also available in a rollator configuration with wheels sold separately.

Dimensions

Length	48cm
Width	75cm
Height	74 - 94cm (Adjustable)
Foot	Fitted with 22mm foot plugs at lower end

Weight Limit

Weight Limit	300kg
--------------	-------

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, such as using appropriate footwear and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.