

Folding Rollator Two Push Button | Technical Spec Sheet

Product Overview

Description: An easy-to-use folding rollator
Colours: Available in silver colour only
Handles: Soft grip handles



Specifications

Description

An easy-to-use folding rollator designed to be used primarily by people who are recuperating from leg or back injuries as well as hip and knee replacements. This rollator can be folded with the push of a quick double release button allowing for efficient storage and travel.

The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy. The two wheels at the front make for easier and more energy efficient propulsion whilst still providing a stable base of support.

Dimensions

Length	58cm
Width	49cm
Height	80 - 87cm (Adjustable)
Wheels	125mm diameter wheels

Weight Limit

Weight Limit	110kg
--------------	-------

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, such as using appropriate footwear and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.